



Cage Washer Testing & Inspection Log: CNTC-001

Date	Facility	Washer #	Rack #	Name

WATER QUALITY: AQUATEST									
Color Chart					Target Values	Cold Water Results	Hot Water Results	Comments	
pH	2.0	4.0	5.0	6.5	8.5	9.5	10.5	12.0	
Alkalinity	0	40	80	120	180	240	500		
Hardness	0	50	120	250	425	1000			

WASHER EFFICIENCY: CAGEWASHCHECK									
Reference	Location	Temperature (Circle)			Water Presence (Circle)		Comments		
 <p>Target Temperature:</p> <p>Place in 4 corners of cage</p> <p>Black ink will wash away if contacted by water</p>	Top Left	120	150	180	Pass	Fail			
	Bottom Left	120	150	180	Pass	Fail			
	Top Right	120	150	180	Pass	Fail			
	Bottom Right	120	150	180	Pass	Fail			

CLEANING EFFICACY: TOSI									
Result Guide									
0									
1									
2									
3									
4									
5									
Result:									
Comments									

PROTEIN DETECTION: PROCHEK-II			
Reference	Pass/Fail	Location Tested of Cage	Comments
 <p>Choose area of cage to swab. A color change to blue-green indicates presence of protein.</p>			

CAGE WASHER INSPECTION		
Inspection	Yes/No	Comments
Spray nozzles are present		
Spray nozzles are free of debris		
Nozzles are aligned, pointed at cage		
Washer screen is free of debris		
Cage washer is free of staining or scaling		
Detergent is at correct levels		
Cycle program is correct		

CART CONDITION	
Area	Comments
Casters	
Shelves	